


How to use the Uniblock


*Versatile support for physiotherapists,
paramedics, and clinical practitioners.*



Step 1 Select the Elevation You Need


Each side of the Uniblock offers a different height or curve. Flip or rotate the block depending on the support required for the procedure.

 **Tip:** The angled sides are ideal for the neck or knees; flat surfaces work best for limbs or general support.



Step 2 Position for Comfort and Stability


Place the Uniblock under the head, arms, knees, legs, or feet for targeted elevation and alignment.

 **Tip:** You can easily flip or reposition the block mid-treatment to adjust to patient movement or changes in positioning needs.



Step 3 Clean After Use

Wipe down with disinfectant spray or clinic wipes between patients.

 **Tip:** For best results, avoid abrasive cleaners and dry thoroughly before storing to maintain hygiene and longevity.



5 Elevations Options

*One block, multiple
heights*



Compact & Portable

Easy to transport



Built to Last

*Polyurethane body with
non-slip, wipeable surface*

Have questions or need support?